

Stellaris Elite

Vision Enhancement System



CATARACT



LASER



RETINA

BAUSCH + LOMB
See better. Live better.

EXPERIENCE THE DUAL-PORT DIFFERENCE

We are constantly seeking new ways to put innovation into your practice, and the Bi-Blade vitrectomy cutter is a perfect example of our dedication to providing precise and efficient surgical experiences.

STABILITY 15,000 CPM

The Bi-Blade cutter's dual-edge blade effectively doubles the cut rate, which allows for faster procedures while potentially decreasing the risk of vitreoretinal traction⁽¹⁾.

EFFICIENCY

Thanks to its unique dual port design that allows for cutting without port closure.

As the port remains open, the Bi-Blade provides consistent aspiration at all times during cutting. The Bi-Blade cutter's have been designed to create a consistent flow rate that operates independently from the cut rate.

1. Teixeira A. et al. Vitreoretinal traction created by conventional cutters during vitrectomy. Ophthalmology 2010;117:1387-1392

Stellaris Elite™ Combined Pack Bi-Blade®

SE5523MVB / SE5523WVB
23G Combined Pack

SE5525MVB / SE5525WVB
25G Combined Pack

SE5527MVB
27G Combined Pack



Stellaris Elite™ Posteriors Pack Bi-Blade®

SE5423MVB / SE5423WVB
23G Posterior Pack

SE5425MVB / SE5425WVB
25G Posterior Pack

SE5427MVB
27G Posterior Pack



Stellaris Elite™ Bi-Blade® Vitrectomy Cutter

SE5623B 23G Vitrectomy cutter
Bi-Blade®

SE5625B 25G Vitrectomy cutter
Bi-Blade®

SE5627B 27G Vitrectomy cutter
Bi-Blade®

